
behind the ear hearing aid, but as I wear glasses I found that the frame and the hearing aid behind my ear constantly clashed and required adjustment. Also the tube from hearing aid to ear canal was uncomfortable on my ear and also very unsightly and my wife said it made me look older. Also, there was an echo. NOW, I wear a red in-ear mini in my right ear and the sound quality is excellent and I hardly know I am wearing a hearing aid.and most people I speak with do not know either. BEST of ALL, my wife has stopped complaining that I have misunderstood what she has said or that she has to repeat herself or shout at me, because I am happy to wear the aid from morning to night. • HOW TO SAVE A FORTUNE: NHS hearing tests are free – if you think you have hearing loss in one or both ears ask your GP to arrange a hearing test. Then, you can accept the behind the ear hearing aids provided by the NHS or if you want in in-ear, less visible aid that will not clash with your glasses or sun-glasses purchase your new in-ear aids from us. The price at Boots, Specsavers, Hidden Hearing and other high-street providers could be £1000 or more. Technical Parameters: • The hearing aid is built to EU standards and bears the CE mark and also is certified that it contains no dangerous materials and has the ROHS mark. • MAX-SSPL90: 118±3 dB SPL • HFA-SSPL90: 110±4 dB SPL • Peak Gain: 33 dB • HFA Full-on Gain: 30±5 dB • Frequency Range: 200-4000 Hz • THD : ?10% • Eq. Input Noise Level: 28 dB SPL • Battery Current Drain: 0.6 mA Each Hearing Aid Package includes: 1 x hearing aid 4 x Different size earplugs 1 x screwdriver 1 x cleaning brush 1 user instruction leaflet 4 packs of 6 A10 Rayovac Hearing Aid